

Tube Type		Common Tests Refer to Austin Pathology Test Directory	Instructions
	Citrate Blue	Coagulation Studies	Citrate tubes must be a full draw Gently invert 4 times to mix
	EDTA Purple	Haematology: FBE, Blood Group G6PD	Gently invert 8 times to mix
	Serum (Gel) Gold	UE, LFT, Mg, Ca, Glucose TFT (x 2 tubes) Hormones Drug Levels, Antibiotics Immunology Tests	Gently invert 8 times to mix
	Serum (Gel) Light sensitive Gold Amber	Bilirubin (SBR)	Gently invert 8 times to mix
	Lithium Heparin No Gel Green	Chromosomes/Karyotype (x2 tubes, minimum 1ml) Galactose Screen	Gently invert 8 times to mix
	Lithium Heparin Gel Mint Green		Gently invert 8 times to mix
	Fluoride Oxalate Grey	Lactate	Gently invert 8 times to mix